


1401ho, 12, Digital-ro 33-gil, Guro-gu, Seoul, 08377, Korea
 www.healcerion.com | Tel: +(82) 70.7582.6326 | E-mail: info@healcerion.com

HC-BC-Temp-VER 1.0
 Printed in Korea

SONON



HEALCERION

Healcerion's Story

Through the use of mobile technology, Healcerion strives to bring affordable, high-quality medical care to clinical and remote locations around the world. Healcerion offers a mobile-based healthcare system that is specifically designed to provide physicians and patients with flexibility and portability. Founded in 2012 by a healthcare provider, Healcerion understands the importance of accurate, timely diagnostic testing results and actionable information. Healcerion is committed to delivering this information to ensure superior patient care and increased satisfaction for both patients and healthcare providers.

Features



Mobility

With our cable-free technology, any urgent situation, including in the emergency room and operating room can easily be accessed while the transducer itself works as a Wi-Fi that makes it possible to perform where there are no internet or stable electric supplies. The light-weight device makes it possible for easy access with simple performance of Wi-Fi features.










User Friendly

In such emergency situations including internal bleedings or fetal checkups, SONON can use any networks, Wi-Fi, 3G, LTE and 5G to share still images (JPEG and DICOM) and video (MPEG4) to other mobile devices between hospitals for assistance.

SONON

Innovative Mobility, Beyond Ultra HD

SONON 500L Preset

-  Thyroid
-  Carotid
-  Breast
-  MSK
-  Vascular
-  Superficial
-  Thorax

Operating mode

- B Mode
- Color Flow Mode
- PW Mode
- M Mode

Applications

- MSK
- US guided procedure
- Early Screening
- F.A.S.T.
- Trauma
- Training

CREATIVE INNOVATION



Linear **300L**



Linear **500L**

*500L image are rendered design and may differ from the actual product.

Linear (300L) Specification	Linear (500L) Specification
<ul style="list-style-type: none"> Frequency : 5&7.5&10 MHz Weight : 370g Soft AP / 2.4 and 5 GHz Depth : 1~ 10 cm Dimensions : 78 x 229 x 38 mm 	<ul style="list-style-type: none"> Frequency : 5 ~12 MHz Weight : 260g (30% ↓) Soft AP / 2.4 and 5 GHz Depth : 1~ 6 cm Dimensions : 57 x 178 x 39 mm (42% ↓)
<ul style="list-style-type: none"> Battery Type : Rechargeable Li-ion Battery Pack Charging Time : 3 Hours (Full) Capacity : 2600 mAh Using Time : 3 Hours (Scan) 12 Hours (Stand-by) 	<ul style="list-style-type: none"> Battery Type : Rechargeable Li-ion Battery Pack Charging Time : Not yet fixed Capacity : 1150 mAh Using Time : 1 Hours (Scan) 8 Hours (Stand-by)

*Product specification above may differ from the actual product.

Item	Minimum requirements
Compatible Device	<ul style="list-style-type: none"> iPhone 6S (2015) or later iPad Air 2 / iPad Mini 4 or later / iPad Pro Galaxy S6, Galaxy Note (2015) or later Galaxy Tab S3 (2017) or later
Platform Free	<ul style="list-style-type: none"> iOS 11.0 or later Windows 10 10.0.15063.0 or later

